



Recombinant *Oryza sativa* subsp. *indica* Homeobox-leucine zipper protein HOX20 (HOX20)

Product Code	CSB-YP386116OFF
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A2YWC0
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>indica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MRSPAALLPV VADGGGGVGV EEEMDVDEDM AMCGGRGGGG GEKKRRLSVE QVRALERSFE TENKLEPERK ARLARDLGLQ PRQVAVWFQN RRARWKTQKL ERDYAALRQS YDALRADHDA LRRDKDALLA EIKELKGKLG DEDAASFSS VKEEEDPAAS DADPPATGAP QGSSES DSSA VLND AEILPH KPAPAAAADA AASEETEAVV TGAALLHHA E VFFHGQLLKV DDDEAAFLGD DGAACGGFFA DEHLPSLPWW AEPTEQWTT
Source	Yeast
Target Names	HOX20
Protein Names	Recommended name: Homeobox-leucine zipper protein HOX20 Alternative name(s): HD-ZIP protein HOX20 Homeodomain transcription factor HOX20 OsHox20
Expression Region	1-269
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.