



Recombinant *Oryza sativa* subsp. *indica* Uncharacterized protein Osl_030282 (Osl_030282)

Product Code	CSB-EP386120OFF-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A2Z139
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>indica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	VLVQ VECDAAVEKS FNKALLAPVD KRLDEAAQAI NEAADSVVAA APPAKKDEVE AATWKRRMFA ITALGMAQGD EKKVAATSLA YKKAACKAVLD AAPADKFKLM DESFKVAVMQ VIAS
Source	<i>E.coli</i>
Target Names	Osl_030282
Protein Names	Recommended name: Uncharacterized protein Osl_030282 Alternative name(s): Unknown protein AN08 from 2D-PAGE of anther
Expression Region	17-134
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.