



Recombinant *Xenopus laevis* Evolutionarily conserved signaling intermediate in Toll pathway, mitochondrial (ecsit)

Product Code	CSB-BP387207XBE
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A2VD95
Product Type	Recombinant Protein
Immunogen Species	<i>Xenopus laevis</i> (African clawed frog)
Purity	>85% (SDS-PAGE)
Sequence	<p> MTCSLHT SQPPGPPQTT SSAGSVAPYE DVFKKEQRNK SNFVEVLDLY CKRDTRRRGH VEIETALRW MPEFGVEKDL EVYNKLLDVF PKEIFIPQNF IQRMFNHYPR QQECAIRVLE QMEHYGITPN KQTCFLLLI FGDRSHPIKK YQRMMYWFPF KYTNPFPVP AELPADPVEL SRMCLHRIAA DLDSRVAVHQ MPYTEKCDDG KEQEHLHIVG IQSPSQVSLI AAHDPSTPVY VEGPFPLWLK KTCVYYYILR ADPDPMENE EIDNERSFYI PLSLDIQLDR DLGDDDSFDV DDVVEGPVFA MCMAGSGDER LLGKWIRGLQ ETNPILGQTP VVFRNLNSGPQ ELSVPQEKEA QPKQEVVEEEE REPERHRMEQ </p>
Source	Baculovirus
Target Names	ecsit
Protein Names	Recommended name: Evolutionarily conserved signaling intermediate in Toll pathway, mitochondrial
Expression Region	34-410
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.