



Recombinant *Oryza sativa* subsp. *indica* Actin-1 (ACT1)

Product Code	CSB-MP387237OFF
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A2XLF2
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>indica</i> (Rice)
Purity	≥85% (SDS-PAGE)
Sequence	MADAEDIQPL VCDNGTGMVK AGFAGDDAPR AVFPSIVGRP RHTGVMVGMG QKDAYVGDEA QSKRGILTLK YPIEHGIVSN WDDMEKIWHH TFYNELRVAP EEHPVLLTEA PLNPKANREK MTQIMFETFN TPAMYVAIQA VLSLYASGRT TGIVLDSGDG VSHTVPIYEG YALPHAILRL DLAGRDLTDY LMKILTERGY SFTTTAEREI VRDMKEKLSY IALDYDQEME TAKTSSSVEK SYELPDGQVI TIGAERFRCP EVLFQPSFIG MEAAGIHETT YNSIMKCDVD IRKDLYGNIV LSGGTTMFPG IADRMSKEIT ALAPSSMKIK VVAPPERKYS VWIGGSILAS LSTFQQMWIA KAEYDESGPS IVHRKCF
Source	Mammalian cell
Target Names	ACT1
Protein Names	Recommended name: Actin-1
Expression Region	1-377
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.