



Recombinant *Oryza sativa* subsp. *indica* Zinc finger A20 and AN1 domain-containing stress-associated protein 8 (SAP8)

Product Code	CSB-EP387259OFF-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A2YEZ6
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>indica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MEHKETGCQQ PEGPILCINN CGFFGSAATM NMCSKCHKEM IMKQEQAKLA ASSIDSIVNG GDSGKEPIIA GHAEVAVAQV EVKTLVAQPA EIAGPSEGVT VNPKGREGPN RCSTCRKRVG LTGFNCRCGN LYCAMHRYSD KHDCQFDYRT AARDAIAKAN PVVKAEKLDK I
Source	E.coli
Target Names	SAP8
Protein Names	Recommended name: Zinc finger A20 and AN1 domain-containing stress-associated protein 8 Short name= OsSAP8 Alternative name(s): Stress-associated protein 3
Expression Region	1-171
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.