



Recombinant *Oryza sativa* subsp. *indica* Homeobox-leucine zipper protein HOX18 (HOX18)

Product Code	CSB-EP387261OFF-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A2YGL9
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>indica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MEGEDLGSWL GLGIGGGGYA YGGDDCRRSP SSPSPVQMLF SQHVKEEITR GYDHGRDEEQ ASGSKIMKGE RGARLRVMRS IRNSGGDGSR SRVLSLGDDG GDGGSGSGGG GGTRKKLQLT KEQSTLLEDS FRVHNILSHA QKHELARQLK LKPRQVEVWF QNRRARTKLK QTEVDCEFLK RCCESLTEEN KQLKHELMEL RRLASAAAAA AGSPLYVQFP RAAAAAMVNV CPSCEKVTVM GGGGGETGKS SSSYSS
Source	E.coli
Target Names	HOX18
Protein Names	Recommended name: Homeobox-leucine zipper protein HOX18 Alternative name(s): HD-ZIP protein HOX18 Homeodomain transcription factor HOX18 OsHox18
Expression Region	1-256
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.