



Recombinant *Oryza sativa* subsp. *indica* Dehydration-responsive element-binding protein 1J (DREB1J)

Product Code	CSB-MP387274OFF
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A2YXQ7
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>indica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MDVARDMEKN TTAMGQLMSS SATTAAATATG PASPKRPAGR TKFQETRHPV FRGVRRRGRA GRWVCEVRVP GSRGDRLWVG TFDTAEEAAR AHDAAMLAMC GASASLNFTD SAWLLHVPRA PVASGHDQLP DVQRAASEAV AEFQRRGSTA ATATATSGDA ASTAPPSSSP VLSPNDDNAS SASTPAVAAA LDHGDMFGGM RADLYYASLA QGLLIEPPPP PTTTEGFCDD EGCGGAEMEL WS
Source	Mammalian cell
Target Names	DREB1J
Protein Names	Recommended name: Dehydration-responsive element-binding protein 1J Short name= Protein DREB1J
Expression Region	1-242
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.