



Recombinant *Oryza sativa* subsp. *indica* Homeobox-leucine zipper protein HOX23 (HOX23)

Product Code	CSB-YP387280OFF
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A2Z734
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>indica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MASNGAAAGA MAPFFPPNFL LQMQQPLPLH HQHLQDHAHG GHGGHLLPP PPPSLSPFLP DLAMDAPPPP MYEASGGDGG GGGAASEDEE DGC GG GGGGG GGEKKRRLSV EQVRTLERSF ESGNKLEPER KAQLARALGL QPRQVAIWVQ NRRARWTKKQ LEKDFDALRR QLDAARAEND ALLSLNSKLH AEIVALKGGGA AAAGGGGSSC RQEAASELIN LNVKETEASC SNRSENSSEI NLDISRPAPP PPPPPANESP VNRGIPFYAS IGRGGAGGVD IDQLLLRGGH SPSPAAVTTP PPPKMELGIT GNGGGADAAA AGAGSFGGLL CGAVDEQPPF WPWADGHHHF H
Source	Yeast
Target Names	HOX23
Protein Names	Recommended name: Homeobox-leucine zipper protein HOX23 Alternative name(s): HD-ZIP protein HOX23 Homeodomain transcription factor HOX23 OsHox23
Expression Region	1-351
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.