



Recombinant Human Zinc finger and BTB domain-containing protein 48 (ZBTB48)

Product Code	CSB-MP026340HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P10074
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MDGSFVQHSV RVLQELNKQR EKGQYCDATL DVGGLVFKAH WSVLACCSHF FQSLYGDGSG GSVVLPAGFA EIFGLLLDFF YTGHLALTSG NRDQVLLAAR ELRVPEAVEL CQSFKPKTSV GQAAGGQSGL GPPASQNVNS HVKEPAGLEE EEVSRTLGLV PRDQEPRGSH SPQRPQLHSP AQSEGPSSLC GKLKQALKPC PLEDKKPEDC KVPPRPLEAE GAQLQGG SNE WEVVVQVEDD GDGDYMSEPE AVLTRRKS NV IRKPCAAEPA LSAGSLAAEP AENRKGTA VP VECPTCHKKF LSKYYLKVHN RKHTGEKPF ECPKCGKCYFR KENLLEHEAR NCMNRSEQVF TCSVCQETFR RRMELRVH MV SHTGEMPYKC SSCSQQFMQK KDLQSHMIKL HGAPKPHACP TCAKCFLSRT ELQLHEAFKH RGEKLFVCEE CGHRASSRNG LQMHIKAKHR NERPHVCEFC SHAFTQKANL NMHLRTH TGE KPFQCHLCGK TFRTQASLDK HNRTH TGERP FSCEFCEQRF TEKGPLL RHV ASRHQEGRPH FCQICGKTFK AVEQLRVHVR RHKGVRKFEC TECGYKFTRQ AHLRRHMEIH DRVENYNPRQ RKLRLNIIED EKMVVVALQP PAELEVGS AE VIVESLAQGG LASQLPGQRL CAEESFTGPG VLEPSLIITA AVPEDCDT
Source	Mammalian cell
Target Names	ZBTB48
Protein Names	Recommended name: Zinc finger and BTB domain-containing protein 48 Alternative name(s): Krueppel-related zinc finger protein 3 Protein HKR3 Zinc finger protein 855
Expression Region	1-688
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.