



Recombinant Danio rerio UPF0575 protein C19orf67 homolog (si:ch211-214c7.5)

| | |
|--------------------------|---|
| Product Code | CSB-EP387631DIL-B |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | A3KNM4 |
| Product Type | Recombinant Protein |
| Immunogen Species | Danio rerio (Zebrafish) (Brachydanio rerio) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MANTEHDELA SLLRETNEAE QSLTDQEAIT ETVEENTTIA SSFGDGTFLW DEDPLSGYTP SCEDTRIMED KISRIEKQLQ YLLSKADEFQ TQLVWSRNCV QDYGFGNIIP MFLETCQPFF SLLESTARNS NPFQPPQSSN IRMQLLQFSQ QLCSRLEQLV LLFASFSCS LEESDPLSVS HFSIGQCQID NMKVSIFRYC CPTPFLASSL TGLYKRMRWN VEREIEGEGD GKTYDSAEFY FLCCEDVIEA ADDQNEDRSS EGESTETDSE SKVSRIWSIG RWVQTPDPN TEDIIDWVLC SLPCGQYKQL LCLGNEEPSS CTATDCLLGA LLSQETHGTF SIIK |
| Source | E.coli |
| Target Names | si:ch211-214c7.5 |
| Protein Names | Recommended name: UPF0575 protein C19orf67 homolog |
| Expression Region | 1-344 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |