



Recombinant Danio rerio UPF0575 protein C19orf67 homolog (si:ch211-214c7.5)

Product Code	CSB-MP387631DIL
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A3KNM4
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MANTEHDELA SLLRETNEAE QSLTDQEAIT ETVEENTTIA SSFGDGTFLW DEDPLSGYTP SCEDTRIMED KISRIEKQLQ YLLSKADEFQ TQLVWSRNCV QDYGFGNIIP MFLETCQPFF SLLESTARNS NPF CPPQSSN IRMQLLQFSQ QLCSRLEQLV LLFASF SFC S LEESDPLSVS HFSIGQCQID NMKVSIFRYC CPTPFLASSL TGLYKRMRWN VEREIEGEGD GKTYDSAEFY FLC CEDVIEA ADDQNE DRSS EGESTETDSE SKVSRIWSIG RWVQTPDPN TEDIIDWVLC SLPCGQYKQL LCLGN EEPSS CTATDCLLGA LLSQETHGTF SIIK
Source	Mammalian cell
Target Names	si:ch211-214c7.5
Protein Names	Recommended name: UPF0575 protein C19orf67 homolog
Expression Region	1-344
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.