



Recombinant Human Myb-related protein A (MYBL1)

Product Code	CSB-MP015263HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P10243
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	<p>MAKRSRSEDE DDDLQYADHD YEVPQQKGLK KLWNRVKWTR DEDKLLKLV EQHGTDDWTL IASHLQNRSD FQCQHRWQKV LNPelikGPW TKEEDQRVIE LVQKYGPKRW SLIAKHLKGR IGKQCRERWH NHLNPEVKKK SWTEEDRII YEAHKRLGNR WAEIAKLLPG RTDNSIKNHW NSTMRRKVEQ EGYLQDGIKS ERSSSKLQHK PCAAMDHMQT QNQFYIPVQI PGYQYVSPEG NCIEHVQPTS AFIQQPFIDE DPDKEKKIKE LEMLLMSAEN EVRRKRIPSQ PGSFSSWSGS FLMDDNMSNT LNSLDEHTSE FYSMDENQPV SAQQNSPTKF LAVEANAVLS SLQTIPEFAE TLELIESDPV AWSDVTSFDI SDAAASPIKS TPVKLMRIQH NEGAMECQFN VSLVLEGKKN TCNGGNSEAV PLTSPNIAKF STPPAILRKK RKMRVGHSPG SELRDGSLND GGNMALKHTP LKTLPFSPSQ FFNTCPGNEQ LNIENPSFTS TPICGQKALI TTPLHKETTP KDQKENVGFR TPTIRRSILG TTPRTPTPFK NALAAQEKKY GPLKIVSQPL AFLEEDIREV LKEETGTDLF LKEEDEPAYK SCKQENTASG KKVRKSLVLD NWEKEESGTQ LLTEDISDMQ SENRFTTSLL MIPLEIHDN RCNLIPEKQD INSTNKTYTL TKKKPNPNTS KVKLEKNLQ SNCEWETVVY GKTEDQLIMT EQARRYLSTY TATSSTSRAL IL</p>
Source	Mammalian cell
Target Names	MYBL1
Protein Names	Recommended name: Myb-related protein A Short name= A-Myb Alternative name(s): Myb-like protein 1
Expression Region	1-752
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself.
Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.