



# Recombinant Rhodobacter sphaeroides

## Fructose-1,6-bisphosphatase class 1 2 (fbp2)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                             |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-MP388542RIJ                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                          |
| <b>Uniprot No.</b>       | A3PRW1                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Immunogen Species</b> | Rhodobacter sphaeroides (strain ATCC 17029 / ATH 2.4.9)                                                                                                                                                                                                                                                                                                                                                     |
| <b>Purity</b>            | ≥85% (SDS-PAGE)                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Sequence</b>          | MAIELEDLGL SPDVADVMDQR LARVGAGIAR IISRNGLERD LGAGVGTNAG<br>GDGQKALDVI ADDAFRAALE GSAVAYYASE EQDEVVTLGE GSLALAIPL<br>DGSSNIDVNV SIGTIFSIFP AAAGPEASFL RPGTEQIAGG YIIYGPQCAL<br>VCSFGQGVQH WVLDLGAGIF RRMPDIRPLP AETSEFAINA SNYRHWPQPI<br>RAFVDDLVAE AEGPRGKNFN MRWIASLVAE THRILMRGGV FLYPGDERKG<br>YERGLRHVY ECAPIAFLIA NVGGRATDGC ADILTALPDR LHARTPFVFG<br>CASKVARVAA YHDLACEETS ALFGSRGLFR S               |
| <b>Source</b>            | Mammalian cell                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Target Names</b>      | fbp2                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Protein Names</b>     | Recommended name: Fructose-1,6-bisphosphatase class 1 2 Short name=FBPase class 1 2 EC= 3.1.3.11 Alternative name(s): D-fructose-1,6-bisphosphate 1-phosphohydrolase class 1 2                                                                                                                                                                                                                              |
| <b>Expression Region</b> | 1-331                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                         |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                                                                               |
| <b>Protein Length</b>    | full length protein                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                          |