



Recombinant Shewanella loihica Sulfate adenylyltransferase subunit 2 (cysD)

Product Code	CSB-MP388641STO
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A3QCU1
Product Type	Recombinant Protein
Immunogen Species	Shewanella loihica (strain ATCC BAA-1088 / PV-4)
Purity	>85% (SDS-PAGE)
Sequence	MSTQQTHLKQ LEAESIQIMR EVAAEFENPV MLYSVGKDSS VLLHLARKAF YPGKIPFPLL HVDTNWKFKE MIEFRDRMAK QYGFDLIVHK NPRGLEMNIS PFTHGSAKHT DIMKTEGLKQ ALDMHGFDAF FGGARRDEEK SRAKERVYSF RDSKHRWDPK NQRPELWNIY NGKVDKGESI RVFPLSNWTE LDIWQYIYLE NIEIPSLYLA EPRPVVKRDG TLIMVDDERM ELEPGEAVEQ KMRVFRFTLGC YPLTGAVESQ ATTLPEIIQE MLLCTTSERQ GRVIDNDSAG SMEKKKMEGY F
Source	Mammalian cell
Target Names	cysD
Protein Names	Recommended name: Sulfate adenylyltransferase subunit 2 EC= 2.7.7.4 Alternative name(s): ATP-sulfurylase small subunit Sulfate adenylate transferase Short name= SAT
Expression Region	1-301
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.