



Recombinant Danio rerio RNA-binding protein 48 (rbm48)

Product Code	CSB-YP389377DIL
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A4FVJ7
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MDGSVASSSV WDTQKVYKHH EQQNV AQTRP KYREGRR LKA VKVY TINLES RYLLVQGVPA IGVMAELVQL FALYGVIEEY RPLDEYP AEQ FTEVYLIQFQ KLTSARA AKR HTDEKSFFGG QLVHVCY APEY ETVEETKQKL QDRRRYV NWA SQNAAKLHSQ QAEVNTESSS STDTRTAEAP IMQKNPEEAR RENVNSDYMG FPLLPLPPTA DISYRLQNHF TQPSQLQWTK ETTE DKMGSL HNFIPPVQKT STQSESSSSS GVKEGDLRQK QKISTPSIRF MPRTT HLESR KRKLDEQ TFF LKEADKTGVL IGPKLPELPK LDMEDSSLNV TANLIRNTMT KAASVPEAKP VQAKHTTPKP RRR I
Source	Yeast
Target Names	rbm48
Protein Names	Recommended name: RNA-binding protein 48
Expression Region	1-364
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.