



Recombinant Magnaporthe oryzae Protein transport protein SEC13 (SEC13)

Product Code	CSB-EP390705MPM
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A4REK3
Product Type	Recombinant Protein
Immunogen Species	Magnaporthe oryzae (strain 70-15 / ATCC MYA-4617 / FGSC 8958) (Rice blast fungus) (Pyricularia oryzae)
Purity	>85% (SDS-PAGE)
Sequence	MIHDAVLDDYY GRRLATCSSD RTIKIFEVEG ETHRLTETLK GHEGAVWCVA WAHPKYGNIL ASSGYDGKVF IWREQGGAWQ KIFDFALHKA SVNIVSWSPH ESGCLLACAS SDGHVSVLEF KDNSFDHQTF LAHGQGVNSV SWAPSTAPGS IISTNATPAA QRRFVTGGSD NTLKIWSWDA ASAQYRCEEG GVLSGHTDWW LDVDWSPTVL QKSYIASASQ DRTVRIWTS DSSNPGVWQSR VLKEFDTTVW RVSWLSLGNV LAVSSGDNKV TLWKENLKGE WACVNSLED
Source	E.coli
Target Names	SEC13
Protein Names	Recommended name: Protein transport protein SEC13
Expression Region	1-289
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.