



# Recombinant *Aeromonas salmonicida* Maltoporin 1 (IamB1)

<b>Product Code</b>	CSB-MP390832AUI
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	A4SNG2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Aeromonas salmonicida</i> (strain A449)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	VDFSG YFRSGVGVVSQ DGDMQTGNQS LVGRLGNESD TYTEIGIGQE VYNKDGKVFY VDSMFSMQSN GSNDYESTAT VCDFDKKQCS EDATFALRQF NVKAKGLISA APDAVVWAGK RFYQRHDLHI IDTKYWNISG AGAGIENLKA GPGAFSLAWI RSDGNDIDNS IVDNDLNVNF LDLRYAGWAP WEGAWTEAGV SYAMPNPTDE QKATGGKYDP ENGVMLTAEM SQYFAGSGIN EKLVLQYANK GLAQNMISQG GGWYDVWQLT DDAKGYRVIL TGDIPLGDKF SVNHVFTYGK GEKLQEWHDN TELFSAVARG GYAWTDIMKT LVEAGTYEST KTWTSGAEDK SSGQKYTLAQ AWSAGPSMFA RPEIRVFASY LKDGEGESFN GGEDDSTWNF GVQAEAWW
<b>Source</b>	Mammalian cell
<b>Target Names</b>	IamB1
<b>Protein Names</b>	Recommended name: Maltoporin 1 Alternative name(s): Maltose-inducible porin 1
<b>Expression Region</b>	26-428
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.