



Recombinant Danio rerio UPF0760 protein C2orf29 homolog (zgc:163002)

Product Code	CSB-BP391619DIL
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A4QP78
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MTLAPKDLSN LLGIMSEDAC SSSFEQLSSA FHHCFGKAEH FGVGSLVVML LQQPDLLPSA PQRLTALYLL WEMYRTEPLA ANPFAAVFAH LLNPSSAAEE QEKTLGFLP PITLPEKFFL SQLMLAPPRE LFKKTPRQVS CMDITSTPQT IDISGLQLAL AERQSELPTQ SKASFPSLLS DPDPDSSNSG FDSSVANQIT ESLVTGPRPP LESHFRPEFI RPTPLLVC E DELCWLNPCE PEHRVQWDRS MCVRNGTGVE IKRIMAKAFK SPLSAAQQTQ LLGELEKDLK LVYHIGLTPA KLPDLVENNP LVAIEMLLKL MSSSQITEYF SVLVNMDMSL HSMEVVNRLT TAVDLPPEFI HLYISNCIST CEQIKDKYMQ NRLVRLVCVF LQSLIRNKII NVQDLFIEVQ AFCIEFSRIR EAAGLFRLK TLDNGETPAE GKTSK
Source	Baculovirus
Target Names	cnot11
Protein Names	Recommended name: UPF0760 protein C2orf29 homolog
Expression Region	1-445
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.