



Recombinant Mouse Retinoic acid receptor alpha (Rara)

| | |
|--------------------------|--|
| Product Code | CSB-YP019338MO |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | P11416 |
| Product Type | Recombinant Protein |
| Immunogen Species | Mus musculus (Mouse) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MASNSSSCPT PGGGHLNGYP VPPYAFFFP MLGGLSPPGA LTSLQHQLPV SGYSTPSPAT IETQSSSSEE IVPSPSPPP LPRIYKPCFV CQDKSSGYHY GVSACEGCKG FFRRSIQKNM VYTCHRDKNC IINKVTRNRC QYCRLQKCFD VGMSKESVRN DRNKKKKEAP KPECSesyTL TPEVGELIEK VRKAHQETFP ALCQLGKYTT NNSSEQRVSL DIDLWDFSE LSTKCIKTV EFAKQLPGFT TLTIADQITL LKAACLDILI LRICTRYTPE QDTMTFSDGL TLNRTQMUNA GFGPLTDLVF AFANQLLPLE MDDAETGLLS AICLICGDRQ DLEQPDKVDM LQEPLLEALK VYVRKRRPSR PHMFPKMLMK ITDLRSISAK GAERVITLKM EIPGSMPLI QEMLENSEGL DTLSGQSGGG TRDGGGLAPP PGSCSPSLSP SSHRSSPATQ SP |
| Source | Yeast |
| Target Names | Rara |
| Protein Names | Recommended name: Retinoic acid receptor alpha Short name= RAR-alpha Alternative name(s): Nuclear receptor subfamily 1 group B member 1 |
| Expression Region | 1-462 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |