



# Recombinant Danio rerio tRNA (uracil (54)-C (5))-methyltransferase homolog-B

|                          |  |
|--------------------------|--|
| <b>Product Code</b>      | CSB-EP393461DIL  |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |
| <b>Uniprot No.</b>       | A4QP75   |
| <b>Product Type</b>      | Recombinant Protein  |
| <b>Immunogen Species</b> | Danio rerio (Zebrafish) (Brachydanio rerio)  |
| <b>Purity</b>            | >85% (SDS-PAGE)  |
| <b>Sequence</b>          | MSSAFKAMRL HLTCPKPKQV INIPQTLFSY LITNQRAINN SCKPVKTSNA<br>KKPKLPPTTKL SWEERLSDTV TPLWRMTYED QLQWKYEHQK KILLKMMKEL<br>SQDPTRAFSD HLNFPLLPIV ASPVRDGYRN KSTFSVNKGI DGNPKTLGFY<br>IGTGKAGNIV CVHADHLLSI PSKHKMVARC YEDFIRLSPL RPCILFDDGG<br>HWREITIRTN STGHTMAIVY FHPQSLTPEE TDIHKAALVE YFTQGPGAIC<br>QLDSLYFQET TMTRCSHEQS QYQLLYGQTH IYEEVLGFKF RISPDSFFQV<br>NREAAEALYK TVAELSQPCV GGTLLDVCCG TGAIGISLSP QMERVIGIEL<br>IEQAVEDAKF NAALNRVCNC EFLAGKAEVV LPDLMGSLSS DGGLTAVVNP<br>SRAGLHYRVV RALRNHSAIR RLVYISCKPD GEAMRNFREL CCSVGLVRRRI<br>TGEAFKPVVA VPVDLFPHTP HCELVLVFER |
| <b>Source</b>            | E.coli   |
| <b>Target Names</b>      | trmt2b   |
| <b>Protein Names</b>     | Recommended name: tRNA (uracil(54)-C(5))-methyltransferase homolog-B EC=2.1.1.35 Alternative name(s): TRM2 homolog B   |
| <b>Expression Region</b> | 1-480  |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.  |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.  |
| <b>Protein Length</b>    | full length protein  |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.  |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |