




# Recombinant Human Annexin A3 (ANXA3), partial

<b>Product Code</b>	CSB-YP001844HU
<b>Relevance</b>	Inhibitor of phospholipase A2, also possesses anti-coagulant properties. Also cleaves the cyclic bond of inositol 1,2-cyclic phosphate to form inositol 1-phosphate
<b>Abbreviation</b>	Recombinant Human ANXA3 protein, partial
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P12429
<b>Alias</b>	35-alpha calcimedlin Annexin III Annexin-3 Inositol 1,2-cyclic phosphate 2-phosphohydrolase Lipocortin III Placental anticoagulant protein III
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥ 90% as determined by SDS-PAGE.
<b>Sequence</b>	WVGHRGTVRDYPDFSPSVDAEAIQKAIRGIGTDEKMLISILTERSNAQRQLIVK EYQAAYGKELKDDLKGDLSGFHEHLMVALVTPPAVFDKQLKKSMMKGAGTNE DALIEILTRTRSRQMKDISQAYYTVYKKS LGDDISSETSGDFRKALLTLADGRRD ESLKVDEHLAKQDAQILYKAGENRWGTDEDKFTEILCLRSFPQLKLTDFDEYRNI SQKDIVDSIKGELSGHFEDLLLAIVNCVRNTPAFLAERLHRALKGIGTDEFTLNRI MVSRSIDLLDIRTEFKKHGYSLYSAIKSDTSGDYEITLLKICGGD
<b>Research Area</b>	Signal Transduction
<b>Source</b>	Yeast
<b>Target Names</b>	ANXA3
<b>Protein Names</b>	Recommended name: Annexin A3 Alternative name(s): 35-alpha calcimedlin Annexin III Annexin-3 Inositol 1,2-cyclic phosphate 2-phosphohydrolase Lipocortin III Placental anticoagulant protein III Short name= PAP-III
<b>Expression Region</b>	5-322aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 6xHis-tagged
<b>Mol. Weight</b>	37.9kDa
<b>Protein Length</b>	Partial
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.