



Recombinant Human Ski-like protein (SKIL)

Product Code	CSB-MP021357HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P12757
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	<p>MENLQTNFSL VQGSTKKLNG MGDDGSPPAK KMITDIHANG KTINKVPTVK KEHLDDYGEA PVETDGEHVK RTCTSVPETL HLNPSLKHTL AQFHLSSQSS LGGPAAFSAR HSQESMSPTV FLPLPSPQVL PGPLLIPSDS STELTQTVLE GESISCFQVG GEKRLCLPQV LNSVLREFTL QQINTVCDEL YIYCSRCTSD QLHILKVLGI LPFNAPSCGL ITLTDAQRLC NALLRPRTFP QNGSVLPKAS SLAQLKETGS AFEVEHECLG KCQGLFAPQF YVQPDAPCIQ CLECCGMFAP QTFVMHSHRS PDKRTCHWGF ESAKWHCYLH VNQKYLGTPE EKCLKIILEE MKEKFMSMRSG KRNQSKTDAP SGMELQSWYP VIKQEGDHVS QTHSFLHPSY YLYMCDKVVA PNVSLTSAVS QSKELTKTEA SKSISRQSEK AHSSGKLQKT VSYPDVSL EE QEKMDLKT SR ELCSRLDASI SNNSTSKRKS ESATCNLVRD INKVGIGLVA AASSPLL VKD VICEDDKGKI MEEVMRTYLK QQEKLNLILQ KKQQQLQMEVK MLSSSKSMKE LTEEQQNLQK ELES LQNEHA QRMEEFYVEQ KDLEKKLEQI MKQKCTCDSN LEKDKEA EYA GQLAELRQRL DHAEADRQEL QDEL RQEREA RQKLEMMIKE LKLQILKSSK TAKE</p>
Source	Mammalian cell
Target Names	SKIL
Protein Names	Recommended name: Ski-like protein Alternative name(s): Ski-related oncogene Ski-related protein
Expression Region	1-684
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.