



# Recombinant Metallosphaera sedula Signal recognition 54 kDa protein (srp54)

<b>Product Code</b>	CSB-MP397423MNM
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	A4YHLO
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Metallosphaera sedula (strain ATCC 51363 / DSM 5348 / JCM 9185 / NBRC 15509 / TH2)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MLDGIRDAVR KFLGGNTSYE VAVEEFIKDI QKALISSDVQ VKLVFSLTNK IKERLKKETP PSNLERREWF IKIVYDELSA LFGGDKEPEV NPKSIPWVIM LVGVQGTGKT TTAGKLAYFY KKRGYKVALV GADVYRPAAL EQLQQIGKQI NVPVYGEPS QDAVGIAKRG VEKFLSERYE LVIVDTAGRH GYGEEVKLLE EMKDIYEKIK PNEVILVIDA SLGQKAYDLA KRFHEASNVG SIIITKMDGT AKGGGALS AV AATGAPIKFI GVGEKIDELE VFNPRRFVAR ILGMDGLEAI IEKMKAMEDY EGIQKKMGEV MTGKSKLTLR DMYKQIVAVR KMGPLSKILQ LIPGMNMMGD IPEDQVKVGE QKMQRWLSIM NSMTYQELDN PSIIDKQRM R RIAMGSGTEV EEVRELIEHF NTVQRTLKML KRRKKDVEKL FGQMGG
<b>Source</b>	Mammalian cell
<b>Target Names</b>	srp54
<b>Protein Names</b>	Recommended name: Signal recognition 54 kDa protein Alternative name(s): SRP54
<b>Expression Region</b>	1-446
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.