



# Recombinant Pan troglodytes Keratin, type I cytoskeletal 17 (KRT17)

<b>Product Code</b>	CSB-MP397684EQV
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	A5A6M0
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Pan troglodytes (Chimpanzee)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MTTIRQFTS SSSIKGSSGL GGGSSRTSCR LSGGLGAGSC RLGSAGGLGS TLGGSSYSSC YSFGSGGGYG SSFSGVDGLL AGGEKATMQN LNDRLASYLD KVRALEEANT ELEVKIRDWY QRQAPGPARD YSQYYRTIEE LQNKILTATV DNANILLQID NARLAADDFR TKFETEQUALR LSVEADINGL RRVLDELTLA RADLEMQIEN LKEELAYLKK NHEEEMNALR GQVGGEINVE MDAAPGVDL RILNEMRDQY EKMAEKNRKD AEDWFFSKTE ELNREVATNS ELVQSGESEI SELRRTMQAL EIELQSLSM KASLEGNLAE TENRYCVQLS QIQGLIGSVE EQLAQLRCEM EQQNQEYKIL LDVKTRLEQE IATYRRLLLEG EDAHLTQYKK EPVTTRQVRT IVEEVQDGKV ISSREQVHQT TR
<b>Source</b>	Mammalian cell
<b>Target Names</b>	KRT17
<b>Protein Names</b>	Recommended name: Keratin, type I cytoskeletal 17 Alternative name(s): Cytokeratin-17 Short name= CK-17 Keratin-17 Short name= K17
<b>Expression Region</b>	1-432
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.