



# Recombinant Danio rerio Iron/zinc purple acid phosphatase-like protein (papl)

<b>Product Code</b>	CSB-BP398839DIL
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	A5D6U8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Danio rerio (Zebrafish) (Brachydanio rerio)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	DVPIGTQ PEQVHISYPG VQNSMLVTWS SANKTDSVVE YGLWGGKLFHS HSATGNSSIF INEGAEYRVM YIHRVLLTDL RPAASYVYHC GSGAGWSELF FFTALNESVF FSPGFALFGD LGNENPQSLS RLQKETQIGT YDVILHIGDF AYDLYEDNGR IGDEFMKQIQ SIAAYVPYMT CPGNHEWAFN FSQYRARFSM PGDTEGLWYS WNVGPAHIIS FSTEYFYFYL EYGLDLLFRQ YEWLRADLQE ANRPENRAER PWIITMGHRP MYCSNDDDDD CTHFQSYVRL GRNDTKPPAP GLEELFYQYG VDLELWAHEH TYERLWPVYD YKVFNGSSEE PYVNPAPVH IITGSAGCRE KHDGFIPKPR DWSAFRSTDY GYTRLQLINN THLYLEQVSD DQYGKVIDQM TLVKEKHGPD AWR
<b>Source</b>	Baculovirus
<b>Target Names</b>	acp7
<b>Protein Names</b>	Recommended name: Iron/zinc purple acid phosphatase-like protein EC=3.1.3.2
<b>Expression Region</b>	24-443
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.