



Recombinant Human Phosphatidylinositol 5-phosphate 4-kinase type-2 beta (PIP4K2B)

Product Code	CSB-MP018023HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P78356
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	SSNCTSTTA VAVAPLSASK TKTKKKHFVC QKVKLFRASE PILSVLMWGV NHTINELSNV PVPVMLMPDD FKAYSKIKVD NHLFNKENLP SRFKFKEYCP MVFRNLRERF GIDDQDYQNS VTRSAPINS SSGRCGTRFL TTYDRRFVIK TVSSEDVAEM HNILKKYHQF IVECHGNTLL PQFLGMYRLT VDGVETVMV TRNVFSHRLT VHRKYDLKGS TVAREASDKE KAKDLPTFKD NDFLNEGQKL HVGEESKKNF LEKLKRDVEF LAQLKIMDYS LLVGIHDVDR AEQEEMEVEE RAEDEECEND GVGGNLLCSY GTPPDSPGNL LSFPRFFGPG EFDPSVDVYA MKSHESSPKK EVYFMAIDI LTPYDTKKKA AHAAKTVKHG AGAEISTVNP EQYSKRFNEF MSNILT
Source	Mammalian cell
Target Names	PIP4K2B
Protein Names	Recommended name: Phosphatidylinositol 5-phosphate 4-kinase type-2 beta EC= 2.7.1.149Alternative name(s): 1-phosphatidylinositol 5-phosphate 4-kinase 2-beta Diphosphoinositide kinase 2-beta Phosphatidylinositol 5-phosphate 4- kinase
Expression Region	2-416
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This protein catalyzes the phosphorylation of phosphatidylinositol-5-phosphate on the fourth hydroxyl of the myo-inositol ring to form phosphatidylinositol-5,4- bisphosphate. This gene is a member of the phosphatidylinositol-5-phosphate 4- kinase family. The encoded protein sequence does not show similarity to other kinases, but the protein does exhibit kinase activity. Additionally, the encoded protein interacts with p55 TNF receptor.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.