



Recombinant Pig Ribonuclease inhibitor (RNH1)

Product Code	CSB-MP019901PI
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P10775
Product Type	Recombinant Protein
Immunogen Species	Sus scrofa (Pig)
Purity	>85% (SDS-PAGE)
Sequence	MNLDIHCEQL SDARWTELLP LLQQYEVVRL DDCGLTEEHC KDIGSALRAN PSLTELCRLT NELGDAGVHL VLQGLQSPTC KIQKLSLQNC SLTEAGCGVL PSTLRSLPTL RELHLSDNPL GDAGLRLLCE GLLDPQCHLE KLQLEYCRLT AASCEPLASV LRATRALKEL TVSNNDIGEA GARVLGQGLA DSACQLETLR LENCGLTPAN CKDLCGIVAS QASLRELDLG SNGLGDAGIA ELCPGLLSPA SRLKTLWLWE CDITASGCRD LCRVLQAKET LKELSLAGNK LGDEGARLLC ESLLQPGCQL ESLWVKSCSL TAACCQHVSL MLTQNKHLL LQLSSNKLGD SGIQELCQAL SQPGTTLRVL CLGDCEVTNS GCSSLASLLL ANRSLRELDL SNNCVGDPGV LQLLGSLEQP GCALEQLVLY DTYWTEEVED RLQALEGSKP GLRVIS
Source	Mammalian cell
Target Names	RNH1
Protein Names	Recommended name: Ribonuclease inhibitor Alternative name(s): Ribonuclease/angiogenin inhibitor 1
Expression Region	1-456
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.