



Recombinant Clostridium botulinum ATP-dependent Clp protease proteolytic subunit (clpP)

Product Code	CSB-MP401627CWV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A5I6W1
Product Type	Recombinant Protein
Immunogen Species	Clostridium botulinum (strain Hall / ATCC 3502 / NCTC 13319 / Type A)
Purity	>85% (SDS-PAGE)
Sequence	MSLVPVVVEQ TNRGERSYDI YSRLKDRII MLSEEVNDTT ASLIVAQLLF LEAEDPKDI HLYINSPGGS ITSGMAIYDT MQYIKPDVST ICVGMAASMG AFLLAAGAKG KRYALPNSEV MIHQPLGGFR GQATDIGIHA ERILKMKKKL NTILSDRTGK PLEQVELDTE RDHFLSAEEA KEYGLIDEVI DKKK
Source	Mammalian cell
Target Names	clpP
Protein Names	Recommended name: ATP-dependent Clp protease proteolytic subunit EC=3.4.21.92 Alternative name(s): Endopeptidase Clp
Expression Region	1-194
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.