



Recombinant Plasmodium vivax Lipoyl synthase, apicoplast (lipA)

| | |
|--------------------------|--|
| Product Code | CSB-MP402469EWW |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | A5K883 |
| Product Type | Recombinant Protein |
| Immunogen Species | Plasmodium vivax (strain Salvador I) |
| Purity | >85% (SDS-PAGE) |
| Sequence | KRVAN AKHATYDMPP KGLRVRDMLE KTAQQNCNQR KRGKCRKFFF LWKMDKMGDT HLGQANGRK NLLRSESATD EASLGGNPLK EKLKESPANW GKDKQEEQQS TDRLPLPKVG NKMPEKKPDW FHVPAPTGKK YNELKADLKK LKLHTVCEEA QCPNIGECWN IGTATIMLLG DTCTRGCKFC SIKTSSKPLP PDANEPFNTA KAICEWDINY VVLTSVDRDD LPDGGASHFA KTIELIKFSR PEILIECLVS DFQGNVDSIR KLANSGMVEVY AHNIETVRRL QKFVRDRRAN YEQSLWVLKT AKEINPLLYT KTSIMLGLGE TKQEVQLQAMA DVRQNNIDVI TFGQYLRPTK NHLNVVEYVS PQMFDYYKEE GMKMGFKYIA SGPLVRSSYK AGEYFMKNLV EQRRGVKLHA EG |
| Source | Mammalian cell |
| Target Names | lipA |
| Protein Names | Recommended name: Lipoyl synthase, apicoplast EC= 2.8.1.8 Alternative name(s): Lipoate synthase Short name= LS Short name= Lip-syn Lipoic acid synthase |
| Expression Region | 26-442 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full Length of Mature Protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |