



Recombinant WU polyomavirus Major capsid protein VP1

Product Code	CSB-BP402662WBG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A5HBD5
Product Type	Recombinant Protein
Immunogen Species	WU polyomavirus (WUPyV)
Purity	>85% (SDS-PAGE)
Sequence	MACTAKPACT AKPGRSPRSQ PTRVQSLPKQ VRKGGVDVLA AVPLSEETEF KVELFVKPVI GNAEGTTPHY WSISSPLKTA EAANVTPDAD TTVCYLSLQV APPDIPNQVS ECDMLIWELY RMETEVLVLP VLNAGILTTG GVGGIAGPQL YFWAVGGQPL DVLGLAPTEK YKGPAQYTVN PKTNGTVPHV YSSSETPRAR VTNEKYSIES WVADPSRNDN CRYFGRMVGG AATPPVVSFS NNSTIPLLDE NGIGILCLQG RLYITCADLL GVNKNRVHTG LSRFFRLHFR QRRVRNPYTI NLLYKQVFNK PADDISGQLQ VTEVTMTEET GPLPPTVEGN VGVPTTSNLS HLPATVTLQA TGPILNTQG
Source	Baculovirus
Protein Names	Recommended name: Major capsid protein VP1 Alternative name(s): Major structural protein VP1
Expression Region	1-369
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.