



Recombinant Danio rerio Otolin-1-A (otol1a)

Product Code	CSB-EP404502DIL-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A5PN28
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	SIDKTTQ WPRMKPTKKP PPRDEGPSKL GSISTTVSPT AIGITEEVD AMMDAYTITS TGSTTFSSDT YSADYHTEAM VPPGVGPGNY TLDYNECFFN FCECCPPERG PPGPVGEKGL PGIPGGKGEM GPPGPPGQEG LTGAPGTHGV KGEKGDGTGAS GLPGIPGVTG KQGEKGESGP KGDKGDGTGFP GLKGDPERG EPGWNGTKGG MGEPGKQGLT GPPGPDGIKG EKGDKGDCPF GEKGQKGSIG EPGPQGPKGD PGVPGTNGTD GLPGSKGPKG DPGPLSKQGE PGPPGPQGPP GQRGMPGMKG TRGLKGARGI RGFKGFKGEP AVQKRSAFSV GLFPSRSFPP PGLPIRFDKI IYNEEAHWDP NASKFNCTHG GYVVFSSYYIT VRNRPLRAL VVNGIRKLRT RDSLYGQDID QASNMAVLRL SSGDQVWLET LRDWNGVYSS SEDDSTFSGF LLYADATKD
Source	E.coli
Target Names	otol1a
Protein Names	Recommended name: Otolin-1-A Alternative name(s): zOtolin1
Expression Region	24-489
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.