



# Recombinant Rabbit Glutathione peroxidase 1 (GPX1)

<b>Product Code</b>	CSB-EP009866RB-B
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P11909
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Oryctolagus cuniculus (Rabbit)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MCAARMAAAA QSVYSFSAHP LAGGEPVNLG SLRGKVLIE NVASLUGTTV RDYTQMNELQ ERLGPRALVV LGFPCNQFGH QENAKNEEIL NSLKYVRPGG GFEPNFMLFQ KCEVNGAKAS PLFAFLREAL PPPSSDPTAL MTDPKFITWC PVCRRNDVSW SFEKFLVGPDPG VPVRRYSRRF PTIDIEPDIQ ALLSKGSGGA
<b>Source</b>	E.coli
<b>Target Names</b>	GPX1
<b>Protein Names</b>	Recommended name: Glutathione peroxidase 1 Short name= GPx-1 Short name= GSHPx-1 EC= 1.11.1.9 Alternative name(s): Cellular glutathione peroxidase
<b>Expression Region</b>	1-200
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Target Details</b>	This gene encodes a member of the glutathione peroxidase family. Glutathione peroxidase functions in the detoxification of hydrogen peroxide, and is one of the most important antioxidant enzymes in humans. This protein is one of only a few proteins known in higher vertebrates to contain selenocysteine, which occurs at the active site of glutathione peroxidase and is coded by UGA, that normally functions as a translation termination codon. In addition, this protein is characterized in a polyalanine sequence polymorphism in the N-terminal region, which includes three alleles with five, six or seven alanine (ALA) repeats in this sequence. The allele with five ALA repeats is significantly associated with breast cancer risk. Two alternatively spliced transcript variants encoding distinct isoforms have been found for this gene.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.