



Recombinant Human Inosine-5'-monophosphate dehydrogenase 2 (IMPDH2)

Product Code	CSB-MP011701HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P12268
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	ADYLISGGT SYVPDDGLTA QQLFNCGDGL TYNDFLILPG YIDFTADQVD LTSALTKKIT LKTPLVSSPM DTVTEAGMAI AMALTGGIGF IHHNCTPEFQ ANEVRKVVKY EQGFITDPVV LSPKDRVRDV FEAARHGFC GIPITDTGRM GSRLVGISS RDIDFLKEEE HDCFLEEIMT KREDLVVAPA GITLKEANEI LQRSKKGKLP IVNEDDELVA IIARTDLKKN RDYPLASKDA KKQLLCGAAI GTHEDDKYRL DLLAQAGVDV VVLDSSQGNS IFQINMIKYI KDKYPNLQVI GGNVVTAQA KNLIDAGVDA LRVGMGSGSI CITQEVLAGC RPQATAVYKV SEYARRFGVP VIADGGIQNV GHIKALALG ASTVMMGSLI AATTEAPGEY FFSDGIRLKK YRGMGSLDAM DKHLSSQNRV FSEADKIKVA QGVSGAVQDK GSIHKFVPYL IAGIQHSCQD IGAKSLTQVR AMMYSGELKF EKRTSSAQVE GGVHSLHSYE KRLF
Source	Mammalian cell
Target Names	IMPDH2
Protein Names	Recommended name: Inosine-5'-monophosphate dehydrogenase 2 Short name= IMP dehydrogenase 2 Short name= IMPD 2 Short name= IMPDH 2 EC= 1.1.1.205 Alternative name(s): IMPDH-II
Expression Region	2-514
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This gene encodes the rate-limiting enzyme in the de novo guanine nucleotide biosynthesis. It is thus involved in maintaining cellular guanine deoxy- and ribonucleotide pools needed for DNA and RNA synthesis. The encoded protein catalyzes the NAD-dependent oxidation of inosine-5 -monophosphate into xanthine-5 -monophosphate, which is then converted into guanosine-5 -monophosphate. This gene is up-regulated in some neoplasms, suggesting it may play a role in malignant transformation.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.