



Recombinant Pseudomonas putida Toluene 1,2-dioxygenase system ferredoxin subunit (todB)

Product Code	CSB-BP405391PZN
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A5W4F0
Product Type	Recombinant Protein
Immunogen Species	Pseudomonas putida (strain ATCC 700007 / DSM 6899 / BCRC 17059 / F1)
Purity	>85% (SDS-PAGE)
Sequence	TWTYILRQG DLPPGEMQRY EGGPEPVMVC NVDGEFFAVQ DTCTHGDWAL SDGYLDGDIV ECTLHFGKFC VRTGKVKALP ACKPIKVFPI KVEGDEVHVD LDNGELK
Source	Baculovirus
Target Names	todB
Protein Names	Recommended name: Toluene 1,2-dioxygenase system ferredoxin subunit
Expression Region	2-107
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.