



Recombinant Human Uncharacterized protein ENSP00000244321

Product Code	CSB-BP405637HU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A6NC86
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	LHCE ICTAAGSRCH GQMKTCSDDK DTCVLLVGKA TSKGKELVHT YKGCIRSQDC YSGVISTTMG PKDHMTSSSF CCQSDGCNSA FLSVPLTNLT ENGLMCPACT ASFRDKCMGP MTHCTGKENH CVSLSGHVQA GIFKPRFAMR GCATESMCFT KPGA EVPTGT NVLFLHHIEC THSP
Source	Baculovirus
Target Names	PINLYP
Protein Names	Recommended name: Uncharacterized protein ENSP00000244321
Expression Region	27-204
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.