



Recombinant Human Ubiquitin carboxyl-terminal hydrolase 17-like protein 3 (USP17L3)

Product Code	CSB-EP406793HU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A6NCW0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MGDDSLYLGGEWQFNHFSKLTSSRPDAAF AEIQR TSLPEKSPLSSETRVDLCD DLAPVAR QLAPREKLPLSSRRPAAVGAGLQNMGNTCYENASLQCLTYTLPLANYMLSRE HSQTCQRP KCCMLCTMQAHITWALHSPGHVIQPSQALASGFHRGKQEDVHEFLMFTVDAM KKACLPGH KQVDHHSKDTTLI HQIFGGCWRSQIKCLHCHGISDTFDPYLDIALDIQAAQSVK QALEQL VKPEELNGENAYHCGLCLQRAPASNTLTLHTSAKVLILVLRFSVDVAGNKLAKN VQYPEC LDMQPYMSQQNTGPLVYVLYAVLVHAGWSCHDGHYFSYVKAQEGQWYKMD DAEVTVC SIT SVLSQQAYVLFYIQKSEWERHSESVSRGREPRALGAEDTDRRAKQGELKRDH PCLQAPEL DEHLVERATQESTLDHWKFLQE QNKTKPEFNVGKVEGTLPPNALVIHQSKYKC GMKNHHP EQQSSLLNLSSTTRTDQESMNTGTLASLQGRTRRAKGKNKHSKRALLVCQ
Source	E.coli
Target Names	USP17L3
Protein Names	Recommended name: Ubiquitin carboxyl-terminal hydrolase 17-like protein 3 EC= 3.4.19.12 Alternative name(s): Deubiquitinating enzyme 17-like protein 3 Ubiquitin thioesterase 17-like protein 3 Ubiquitin-specific-processing protease 1
Expression Region	1-530
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.