



# Recombinant *Klebsiella pneumoniae* subsp. *pneumoniae* Maltoporin 1 (IamB1)

<b>Product Code</b>	CSB-EP407578KAX-B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	A6T5L4
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Klebsiella pneumoniae</i> subsp. <i>pneumoniae</i> (strain ATCC 700721 / MGH 78578)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	ATAAIPT IDFHGYMRAG VGVSGDGSEA EWQKNKLGRL GNESDITYGEL ELGSEVYKGD DVSFYLD SMVSDGSNDN ETTLNDDAQF GLRQLNLQIK GLIPGDPNAV IWGGKRYR HDLHIIDTKY WNISGSGAGV ENYTLGPGAV SLAWIRGDAN DVDYRVDGDS NVNINYIDLR YAGWKPWAGS WTEFGIDYAM PNTTKKQDSY GGLYDADNGV MLTGEISQDM LGGYNKTVLQ YANKGLAQNM VSQGGGWYDM WNYVNDATGY RVINTGLIPI TEKFSINHVL TWGSADDITD YTDKTRMLSL VARGQYQFTD YVRLIGEVGG FYQKDSYNG TSYKQAGEKY TIALGLADGP DFMSRPELRI FASYLNDS ED GKPFE DQTAN NTWNFGVQVE AWW
<b>Source</b>	E.coli
<b>Target Names</b>	IamB1
<b>Protein Names</b>	Recommended name: Maltoporin 1 Alternative name(s): Maltose-inducible porin 1
<b>Expression Region</b>	24-423
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.