



Recombinant *Saccharomyces cerevisiae* ATP-dependent RNA helicase FAL1 (FAL1)

Product Code	CSB-EP408164STA-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A6ZXY5
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain YJM789) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MSFDREEDQK LKFKTSKCLK VSSTFESMNL KDDLRLGIYS YGFEAPSSIQ SRAITQIIISG KDVIAQAQSG TGKTATFTIG LLQAIDLRKK DLQALILSPT RELASQIGQV VKNLGDYMNV NAFITGGKT LKDDLKMKQK HGCQAVSGTP GRVLDMIKKQ MLQTRNVQML VLDEADELLS ETLGFKQQIY DIFAKLPKNC QVVVVSATMN KDILEVTRKF MNDPVKILVK RDEISLEGIK QYVVNVKKEE WKFDTLCDIY DSLTITQCVI FCNTKKKVDW LSQRLIQSNF AVVSMHGDMK QEERDKVMND FRTGHSRVL ISTDVWARGID VQQVSLVINY DLPEIENYI HRIGRSGRFG RKGVAINFIT KADLAKLREI EKFYSIKINP MPANFAELS
Source	E.coli
Target Names	FAL1
Protein Names	Recommended name: ATP-dependent RNA helicase FAL1 EC= 3.6.4.13 Alternative name(s): Translation initiation factor four A-like protein 1
Expression Region	1-399
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.