



Recombinant Human Protein FAM196B (FAM196B)

Product Code	CSB-MP408253HU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A6NMK8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAQQNMKVRP VLLKRNSLES VEFVKQPHHR RSKSQQVRFK EDGTTKNPTG LAEVDVQTPE DPAVMGKTQA TRHHLPTYS LSFPRSQKAG GFRNIAIQTS PSLRKHFPVF KRKRLTASKS LVEMPTASQS AIQVNGNLSE QDIVSSDLAY LRLAQHLEDG PRRVKVSHAF LPRVPKVQSN GPVSICLEAG TWRSLEKATA AIQVPDDIYH SPSWEARESA LSPDRSAEVS NSIHPLDDTR PGDGRRVTPL DSEKSTSCLN ATSVASHTPG TEELKPELLL PKDNSDDKDL GSLSSQSKET CVPSSPRTHS SPSQGSQSHSQP AHPGRASDCP SSSNNHQNLV SLKTNASASKS APGCQEQTAN NPTESDTLEF PNCPGSNHLP SSLSRSETKL QSNREISDIN QIHLARGEIC DLQGRLQSVE ESLHSNQEKI KVLLNVIQDL EKARALTEGR NFYRTGQDLN NCSTCQNTAC IISVEYDFR QQEGRFHEVL QSLEEAEPVE EASPPPKSPA EPPAPEKQDL RRTKKVKKKK CFWWI
Source	Mammalian cell
Target Names	FAM196B
Protein Names	Recommended name: Protein FAM196B
Expression Region	1-535
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.