



Recombinant Human Protein disulfide-isomerase A4 (PDIA4)

Product Code	CSB-YP017722HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P13667
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	VAGAEGPDED SSNRENAIED EEEEEEDDD EEEDDLEVKE ENGLVLNDA NFDNFVADKD TVLLEFYAPW CGHCKQFAPE YEKIANILKD KDPPIPVAKI DATSASVLAS RFDVSGYPTI KILKKGQAVD YEGSRTQEEI VAKVREVSQP DWTTPPEVTL VLTKENFDEV VNDADIILVE FYAPWCGHCK KLAPEYEKAA KELSKRSPPI PLAKVDATEA TDLAKRFDVS GYPTLKIFRK GRPYDYNGPR EKYGIVDYMI EQSGPPSKEI LTLKQVQEFL KDGDDVIIIIG VFKGESDPAY QQYQDAANNL REDYKFHHTF STEIAKFLKV SQGQLVVMQP EKFAQSKYEPR SHMMDVQGST QDSAIFDFVL KYALPLVGHR KVSNDKRYT RRPLVVVYYS VDFSFYRAA TQFWRSKYLE VAKDFPEYTF AIADEEDYAG EVKDLGLSES GEDVNAAILD ESGKKFAMEP EEFSDTLRE FVTAFKKGKL KPVIKSQPVP KNNKGPVKVV VGKTFDSIVM DPKKDVLIIEF YAPWCGHCKQ LEPVYNSLAK KYKGQKGLVI AKMDATANDV PSDRYKVEGF PTIYFAPSGD KKNPVKFEGG DRDLEHLSKF IEEHATKLSR TKEEL
Source	Yeast
Target Names	PDIA4
Protein Names	Recommended name: Protein disulfide-isomerase A4 EC= 5.3.4.1 Alternative name(s): Endoplasmic reticulum resident protein 70 Short name= ER protein 70 Short name= ERp70 Endoplasmic reticulum resident protein 72 Short na
Expression Region	21-645
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.