



# Recombinant *Saccharomyces cerevisiae* Ribosomal N-lysine methyltransferase 5 (RKM5)

<b>Product Code</b>	CSB-YP409779STA
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	A7A136
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Saccharomyces cerevisiae</i> (strain YJM789) (Baker's yeast)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MAFKLWLLDE ETIYEHVFER YTQLEGQSGK LAQDLGIQDR RGGVLEITFE PSGLEGGRKK KRVRRRNKAS SVEEDQNVAV DSYHVSVGQS ISSLRSSRDN GNSTTGYVLW STTPFFINWL LYSTSAAPFR LGSQVEVTCG SSCEGHKLEL PRLVDLTGAD RGKRGILELG AGISGILPVI LGNFVDITYVS TDQKGILNKL KDNIMENLSQ LTRKRCISRS LRLELPTVEP VGDADITAAS LPSKSTLHLE VAALDWEKIN LQDKKTHSLH PELSLIGETC SSVYVIAMDV IYNEYLIDPF LKTLKQLKHW LQTTYNLQFH VLVGIHLRSQ EVTTLFLEKA IIEYDFTVYD IVDQVIQESR FNFYLIT
<b>Source</b>	Yeast
<b>Target Names</b>	RKM5
<b>Protein Names</b>	Recommended name: Ribosomal N-lysine methyltransferase 5 EC= 2.1.1.-
<b>Expression Region</b>	1-367
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.