



Recombinant Human Annexin A8 (ANXA8)

Product Code	CSB-EP001849HU-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P13928
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAWWKSWIEQ EGVTVKSSSH FNPDPDAETL YKAMKGIGTN EQAIIDVLTK RSNTQRQQIA KSFKAQFGKD LTETLKSLS GKFERLIVAL MYPPYRYEAK ELHDAMKGLG TKEGVIIIEIL ASRTKNQLRE IMKAYEEDYG SSLEEDIQAD TSGYLERILV CLLQGSRDDV SSFVDPGLAL QDAQDLYAAG EKIRGTDEM FITILCTRSA THLLRVFEEY EKIANKSIED SIKSETHGSL EEAMLTVVKC TQNLHSYFAE RLYYAMKGAG TRDGLTIRNI VSRSEIDLNL IKCHFCKMYG KTLSSMIMED TSGDYKNALL SLVGSDP
Source	E.coli
Target Names	ANXA8
Protein Names	Recommended name: Annexin A8 Alternative name(s): Annexin VIII Annexin-8 Vascular anticoagulant-beta Short name= VAC-beta
Expression Region	1-327
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Target Details	This gene encodes a member of the annexin family of evolutionarily conserved Ca ²⁺ and phospholipid binding proteins. The encoded protein may function as an anticoagulant that indirectly inhibits the thromboplastin-specific complex. Overexpression of this gene has been associated with acute myelocytic leukemia. A highly similar duplicated copy of this gene is found in close proximity on the long arm of chromosome 10.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.