



# Recombinant Human Annexin A8 (ANXA8)

<b>Product Code</b>	CSB-EP001849HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P13928
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MAWWKSWIEQ EGVTVKSSSH FNPDPDAETL YKAMKGIGTN EQAIIDVLTK RSNTQRQQIA KSFKAQFGKD LTETLKSELS GKFERLIVAL MYPPYRYEAK ELHDAMKGLG TKEGVIIIEIL ASRTKNQLRE IMKAYEEDYG SSLEEDIQAD TSGYLERILV CLLQGSRDDV SSFVDPGLAL QDAQDLYAAG EKIRGTDEMK FITILCTRSA THLLRVFEEY EKIANKSIED SIKSETHGSL EEAMLTVVKC TQNLHSYFAE RLYYAMKGAG TRDGLTIRNI VSRSEIDLNL IKCHFCKMYG KTLSSMIMED TSGDYKNALL SLVGSDP
<b>Source</b>	E.coli
<b>Target Names</b>	ANXA8
<b>Protein Names</b>	Recommended name: Annexin A8 Alternative name(s): Annexin VIII Annexin-8 Vascular anticoagulant-beta Short name= VAC-beta
<b>Expression Region</b>	1-327
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Target Details</b>	This gene encodes a member of the annexin family of evolutionarily conserved Ca <sup>2+</sup> and phospholipid binding proteins. The encoded protein may function as an anticoagulant that indirectly inhibits the thromboplastin-specific complex. Overexpression of this gene has been associated with acute myelocytic leukemia. A highly similar duplicated copy of this gene is found in close proximity on the long arm of chromosome 10.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.