



# Recombinant *Saccharomyces cerevisiae* Autophagy-related protein 14 (ATG14)

<b>Product Code</b>	CSB-BP411324STA
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	A6ZL72
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Saccharomyces cerevisiae</i> (strain YJM789) (Baker's yeast)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MHCPICHHRA HVVYCAHCIN TSPSLLLK LKLDLILLKDEN KELNGKVEQI LNEAMNYDQL DIKRMEKKKD PLMNSLMKLD VLRMCKNNNL IRHRIEQLNE RIYSKRNHIS ELKVEIDNYK CYKVGTTGTDK LIEQVEISDA KNKLAQVSKI CESVRDYKLN LLNNWFVIQK LQDNFQIPFA IAFQPLISLK NFRVLPLAIT NDSINIMWKY ISFFSDILMI KLPYTNKICE QPMFEFSDSI QTVVQRLIKL IINILQICRH LKLVPTPMD IPWLLDQYDV DGLFYNMVKR NKMKCRSVSL YWTFGMLYSM VLDNMNPNQR GHPARRTAPP PTVTGP HDRW YVVG
<b>Source</b>	Baculovirus
<b>Target Names</b>	ATG14
<b>Protein Names</b>	Recommended name: Autophagy-related protein 14 Alternative name(s): Cytoplasm to vacuole targeting protein 12
<b>Expression Region</b>	1-344
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.