



Recombinant *Saccharomyces cerevisiae* Putative cyclin-dependent kinase inhibitor SPL2 (SPL2)

Product Code	CSB-EP411359STA-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A6ZT44
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain YJM789) (Baker's yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MRTYTPLIYN IYNVHIWVFT ESQQQIGQMS PRGKMETAVS QGQHKQLKDG HQHKGRKLSE EIASLLRLKE CRRLNPASY TPRRTSQSQS LSGSTFKEYN EYVNEKDSSR AQRQNAAAVL SKLAHDFWEN DCVIDEDIFE DSSDEEQS
Source	E.coli
Target Names	SPL2
Protein Names	Recommended name: Putative cyclin-dependent kinase inhibitor SPL2 Alternative name(s): Suppressor of PLC1 deletion protein 2
Expression Region	1-148
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.