



Recombinant Human Enolase-like protein ENO4 (ENO4)

Product Code	CSB-EP411467HU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A6NNW6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MEEEGGGRSC GTTRELQKLIK QQAMEYYREN DVPRRLEELL NSTFYLQPAD VYGHKLANCF SKLAKPPTIC KIVGKDVLDG LGLPTLQVDI FCTIQNFPKN VCSVVISTHF EVHENALPEL AKAEEAERAS AVSTAVQWVN STITHELQGM APSDQAEVDH LLRIFFASKV QEDKGRKELE KSLEYSTVPT PLPPVPPPPP PPPPTKKKGQ KPGRKDTITE KPIAPAEPVE PVLSGSMAIG AVSLAVAKAC AMLLNKPLYL NIALLKHNQE QPTTLSMPLL MVSLVSCGKS SSGKLNLMKE VICIPHELT TKQGVEMLME MQKHINKIIE MMPPSPPKAE TKKGHDGSKR GQQQITGKMS HLGCLTINCD SIEQPLLLIQ EICANLGLEL GTNLHLAINC AGHELMODYNK GKYEVIMGTY KNAAEEMVDLY VDLINKYPSI IALIDPFRKE DSEQWDSIYH ALGSRCYIIA GTASKSISKL LEQGNISIPK SNGLIKHHTN QTTMSDLVEI TNLIDSKKHI TVFGSTEGES SDDSLVDLAV GLGVRFIKLG GLSRGERVTK YNRLLTIEEE LVQNGTLGFK EEHTFFYFNE EAEKAAEALE AAAAREPLVP TFPTQGVEES AETGASSG
Source	E.coli
Target Names	ENO4
Protein Names	Recommended name: Enolase-like protein ENO4
Expression Region	1-628
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.