



# Recombinant Human Enolase-like protein ENO4 (ENO4)

<b>Product Code</b>	CSB-EP411467HU
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	A6NNW6
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MEEEGGGRSC GTTRELQKLIK QQAMEYYREN DVPRRLEELL NSTFYLQPAD VYGHKNCNF SKLAKPPTIC KIVGKDVLGD LGLPTLQVDI FCTIQNFPKN VCSVVISTHF EVHENALPEL AKAEEAERAS AVSTAVQWVN STITHELQGM APSDQAEVDH LLRIFFASKV QEDKGRKELE KSLEYSTVPT PLPPVPPPPP PPPPTKKKGQ KPGRKDTITE KPIAPAEPVE PVLSGSMAIG AVSLAVAKAC AMLLNKPLYL NIALLKHNQE QPTTLSMPLL MVSLVSCGKS SSGKLNLMKE VICIPHELT TKQGVEMLME MQKHINKIIE MMPPSPPKAE TKKGHDGSKR GQQQITGKMS HLGCLTINCD SIEQPLLIQ EICANLLEL GTNLHLAINC AGHELMDYNK GKYEVIIMGTY KNAEMVDLY VDLINKYPSI IALIDPFRKE DSEQWDSIYH ALGSRCYIIA GTASKSISKL LEQGNISIPK SNGLIKHHTN QTTMSDLVEI TNLIDSKKHI TVFGSTEGES SDDSLVDLAV GLGVRFIKLG GLSRGERVTK YNRLLTIEEE LVQNGTLGFK EEHTFFYFNE EAEKAAEAL AAAAREPLVP TFPTQGVVEES AETGASSG
<b>Source</b>	E.coli
<b>Target Names</b>	ENO4
<b>Protein Names</b>	Recommended name: Enolase-like protein ENO4
<b>Expression Region</b>	1-628
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.