




# Recombinant Rat Plasma kallikrein (Klkb1), Partial

<b>Product Code</b>	CSB-EP012461RA
<b>Relevance</b>	The enzyme cleaves Lys-Arg and Arg-Ser bonds. It activates, in a reciprocal reaction, factor XII after its binding to a negatively charged surface. It also releases bradykinin from HMW kininogen and may also play a role in the renin-angiotensin syst by converting prorenin into renin.
<b>Abbreviation</b>	Recombinant Rat Klkb1 protein, partial
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P14272
<b>Alias</b>	Fletcher factorKininogenin;Plasma prekallikrein
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Rattus norvegicus (Rat)
<b>Purity</b>	≥ 90% as determined by SDS-PAGE.
<b>Sequence</b>	IVGGTNSSLGEWPWQVSLQVKLVSQNHMCGGSIIGRQWILTAHCFDGIYPYD VWRIYGGILNLSEITNKTPFSSIKELIIHQYKMSEGSYDIALIKLQTPLNYTEFQK PICLPSKADTNTIYTNCWVTGWGYTKERGETQNILQKATIPLVPNEECQKKYRD YVITKQMICAGYKEGGIDACKGDSGGPLVCKHSGRWQLVGITSWGEGCARKE QPGVYTKVAEYIDWILEKIQSSKERALETSPA
<b>Source</b>	E.coli
<b>Target Names</b>	Klkb1
<b>Expression Region</b>	391-638aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 6xHis-tagged
<b>Mol. Weight</b>	31.8kDa
<b>Protein Length</b>	Plasma kallikrein light chain

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.