



Recombinant *Saccharomyces cerevisiae* Increasing suppression factor 1 (ISF1)

Product Code	CSB-EP412204STA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A6ZME0
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain YJM789) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MIASEIFERG VQDPFCQDCD YEDETDVQSF LGSNDLNDV NSKLASFSFE NSSKSNNSHH SSSTNAGNTS RHIGNHTIGH HLRKIKTAPH HLYGFVPASS TNNSNEPIRP SPRRIRANSS TLIHQLSRQS TRQSSLGDAA DSCFDHKCIK PRSRHSSCYG IPHLYGLEK YVSSELDLA VANDQSNDLT SPLTSVSTPA SNSNSYLNLN SSSAAYPSSH LSNEKNNRLK LISHGKISSN NVPGHSGNLN HYHRERTPSN LRRESFLLS NGSSSSPLQT RNNYSNSLV KSPSNSSLNT SVASSNEESI PHTSNCLEER NPRRKSFIKL SLASSFSN
Source	E.coli
Target Names	ISF1
Protein Names	Recommended name: Increasing suppression factor 1 Alternative name(s): Mitochondrial biogenesis regulation protein 3
Expression Region	1-338
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.