



Recombinant *Saccharomyces cerevisiae* Mitochondrial protein import protein ZIM17 (ZIM17)

Product Code	CSB-MP412234STA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A6ZSH0
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain YJM789) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	NEVKKDDKKVHLGSEFKVDKPKMMIAFTCKKCNTRSSHTMSKQAYEKGTVLISC PHCKVRH LIADHLKIFHDHHVTVEQLMKANGEQVSQDVGDLFEFEDIPDSLKDVLGKYAKN NSENASQ LPHPSQK
Source	Mammalian cell
Target Names	ZIM17
Protein Names	Recommended name: Mitochondrial protein import protein ZIM17 Alternative name(s): Mitochondrial import inner membrane translocase subunit TIM15 mtHsp70 escort protein 1 mtHsp70-associated motor and chaperone protein TIM15/ZIM17 Short n
Expression Region	48-174
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.